



Biospin – 유전체R
2015.6.18



Chapter 2

DNA Techniques

Sejin Park





1. DNA 분리기술

- 1) DNA 분리용 시료의 선택
- 2) 항응고제의 선택
- 3) DNA분리의 기본원리
- 4) DNA정량 및 정도 관리
- 5) DNA 시료의 전기영동 분리
- 6) Souther blot 분석

2. 재조합 DNA 클로닝 기술

- 1) DNA조작효소
 - ① DNA복제 효소
 - ② Nuclease
 - ③ DNA결합효소
 - ④ DNA 말단 조작효소
 - ⑤ 역전사효소
- 2) DNA클로닝벡터
- 3) DNA클로닝 기술





3. PCR 증폭기술

- 1) PCR의 원리
- 2) Reverse Transcriptase(RT)-PCR 기술
- 3) Real-time PCR 기술
- 4) Digital PCR 기술





1. DNA 분리기술






DNA 분리용 시료의 선택



- ▶ genomic DNA 분리용 세포 - 조직, 모근, 구강상피세포, 배양된 세포 및 혈액.
- ▶ 많은 양의 DNA 필요 시
 - ▶ 사람의 혈액에 존재하는 B세포를 Epstein-Barr virus(EBV)로 감염시켜 불멸화를 유도 lymphoblastoid cell line(LCLs)을 만들어 사용.
 - ▶ 단점: 세포주 만드는 과정 또는 세포 배양 과정에서 일부 세포주에서 유전체 변이가 발생하는 문제점.
 - ▶ 하지만 대규모 유전체 연구에 사용된 시료는 대부분 EBV-transformed B cell 사용.
- ▶ 하지만!!! 최근에는
 - ▶ whole-genome amplification기술의 발전과 array-based DNA genotyping 분석 기술의 발전으로 적은 양의 유전체 시료로도 많은 종류의 유전체 분석이 가능하게 되었다.
 - ▶ 사람의 혈액으로 부터 바로 DNA를 분리하는 방법을 많이 사용.



Workflow: DNA Preparation to Applications





항응고제의 선택



▶ EDTA

- ▶ DNA분리용 혈액 채취시 사용하는 항응고제
- ▶ DNase 저해기능

▶ heparin

- ▶ DNA분리과정에서 DNA와 결합
 - ▶ PCR증폭에 사용되는 Taq pol을 저해
- ▶ 상온 수송, 채취한지 3일 이상 되면 DNA가 분해된다



○○○ DNA분리의 기본 원리 ○○○

▶ DNA분리과정

- ▶ 세포 파쇄(Cell lysis)
- ▶ lysis
- ▶ 단백질 제거
- ▶ DNA 수집


▶ from 혈액

- ▶ 적혈구 세포 제거 (핵x) 백혈구 세포만을 분리
- ▶ 분리된 백혈구 세포 - 세포막을 파괴
- ▶ 단백질 제거
- ▶ 용액속에 존재하는 DNA를 에탄올로 침전

▶ 다양한 KIT

- ▶ cell lysis : 세포를 lysis buffer로 파괴.
- ▶ 컬럼을 통과 시켜 컬럼에 DNA결합시킨다.
- ▶ 다른 물질들을 씻어낸다.
- ▶ 컬럼에 붙어있던 DNA를 TE buffer로 씻어내어 순수한 DNA만을 분리
- ▶ 분리된 DNA -70도씨에 보관


DNA Isolation Using Commercial Kit



○○○ DNA 정량 및 정도 관리 ○○○

▶ DNA 정량(quantification)

- ▶ UV spectrophotometer
- ▶ nano drop
- ▶ Qubit



▶ $1 A_{260} = 50 \text{ ug/ml}$, $A_{260} = 0.1 \sim 1.0$

▶ A_{260}/A_{280}


= 1.8~2.0 ; 양질의 DNA

< 1.6 단백질 오염 (재침전 필요)

>2.0 RNA오염



DNA Quantitation
Using a Spectrophotometer
<https://www.youtube.com/watch?v=qw2ZaUXgWHU>



Qubit (invitrogen)

qubit 1.0


<https://www.youtube.com/watch?v=RRKZN-7jqg>

qubit 2.0

<https://www.youtube.com/watch?v=bSSIO2fqEN8>


qubit 3.0

<https://www.youtube.com/watch?v=ymny7xABaa4>



nanodrop

<https://www.youtube.com/watch?v=FiGZnNs2xXY>



DNA Quantity and Quality

- **Measured by UV spectroscopy:**



- NanoDrop ND-1000 spectrophotometer (NanoDrop technologies, USA)
 - *particularly good in case of the limited amount of DNA available
 - *it requires 1 ul of samples for measurement

- **Quantification of DNA amount:**

1 A_{260} = 50 $\mu\text{g/ml}$ DNA (range of A_{260} for accurate DNA quantification= 0.1~1.0)

- 260/280 ratio = 1.8 ~ 2.0 → good quality DNA
- 260/280 ratio < 1.6 → protein contamination (should be re-precipitated)
- 260/280 ratio > 2.0 → RNA contamination

Agarose Gel Electrophoresis



Agarose concentration:

high → good for separation of small DNA fragment

low → good for separation of large DNA fragment


→ slow migration

→ fast migration




전기영동



1. Make gel.




2. Obtain prepared DNA samples.




3. Load samples into gel.




4. Separate fragments by electrophoresis.




6. Prepare a standard curve.
Determine fragment sizes.




5. Stain DNA fragments and measure distances.





- 1 Load samples of DNA fragments into wells at the top of the gel.




- 2 Apply an electric field.




- 3 Wait additional time.







Agarose gel을 이용한 DNA의 전기영동




Gel electrophoresis



EtBr Stain & UV light

Southern Hybridization/Blotting




Southern blot	DNA detected by labeled DNA
Northern blot	mRNA detected by labeled DNA
Western blot	Protein detected by antibody (Ab)



2. 재조합 DNA 클로닝 기술



Recombinant DNA Technology





DNA조작효소





1. DNA복제 효소
2. Nuclease
3. DNA결합효소
4. DNA 말단 조작효소
5. 역전사효소



1. DNA복제 효소


DNA- and Primer-dependent DNA Synthesis




- DNA synthesis requires a primer (No DNA synthesis without primer)
- The primers determines which part of a DNA molecule is copied

Multifunctional Activity of DNA Polymerase

(A) 5' → 3' DNA synthesis




(B) 3' → 5' exonuclease activity



proofreading activity
(*pfu enzyme*)

(C) 5' → 3' exonuclease activity



Less common
(*in some DNA polymerases*)

DNA Polymerases

Types of DNA Polymerase	Description	Main Use
DNA polymerase I	Unmodified <i>E. coli</i> DNA polymerase	DNA labeling
<u>Klenow</u> polymerase	Modified <i>E. coli</i> DNA polymerase I	DNA labeling
<u>Sequenase</u>	Modified phage T7 DNA polymerase I	DNA sequencing
<u>Taq</u> polymerase	<u>Thermus aquaticus</u> DNA polymerase I	PCR
Reverse transcriptase	RNA-dependent DNA polymerase (from retroviruses)	<u>cDNA</u> synthesis

E. coli DNA polymerase= polymerase activity
 5' → 3' exonuclease activity
 3' → 5' exonuclease activity

Klenow polymerase= polymerase activity
 3' → 5' exonuclease activity
Lack of 5' → 3' exonuclease activity: good for sequencing/labeling

Optimum Rx temperature for all DNA polymerase except Taq: 37 °C (72 °C for Taq)
 RTase (reverse transcriptase): RNA-dependent DNA polymerase (retroviruses such as HIV)

2. Nuclease

Nucleases

Nuclease	Description	Main Use
Restriction <u>endonucleases</u>	Sequence-specific DNA <u>Endonucleases</u>	Many applications
Nuclease S1	<u>Endonuclease</u> specific for single-stranded DNA and RNA	Transcript mapping
<u>Deoxyribonuclease 1</u>	<u>Endonuclease</u> specific for double-stranded DNA and RNA	Nuclease <u>footprinting</u>


Restriction Enzymes

Enzyme	Recognition sequence	Type of ends	End sequence
<u>AluI</u>	5'-AGCT-3' 3'-TCGA-5'	Blunt end	5'-AG CT-3' 3'-TC GA-5'
<u>EcoRI</u>	5'-GAATTC-3' 3'-CTTAAG-5'	Sticky end (5' overhang)	5'-G AATTC-3' 3'-CTTAA G-5'
<u>PstI</u>	5'-CTGCAG-3' 3'-GACGTC-5'	Sticky end (3' overhang)	5'-CTGCA G-3' 3'-G ACGTC-5'

- >2500 type II enzymes have been isolated
- >300 are available for use in the lab.

Most recognition sequence have inverted symmetry (= *Palindromes*); when read in the 5' → 3' direction, the sequence is the same in both strands.

DNA Ligase



- 3. DNA 결합효소
- 4. DNA 말단 조작 효소
- 역전사효소

Cloning Vectors


1. Plasmid vector
2. Lambda phage vector
3. Bacteriophage vector
4. BACs (bacterial artificial chromosome)
5. YACs (yeast artificial chromosome)
- +
6. Expression vector (for cDNA expression cloning)

1~5: depending on insert size

Bacterial Cloning Systems


Vector type	Cloning capacity (kb)	Examples	Year of appearance	Mode of introduction to E. coli host cell	Mode of growth on agar plate	No. of copies per cell
Plasmids	0.1-12	pBR322, pUC	1977	transformation/ electroporation	colonies	50-100
M13	0.1-3.5	M13mp18	1981	transfection	plaques	N/A
Bacteriophage λ	10-20	EMBL, FIX	1983	packaging and infection	plaques	N/A
Cosmid	35-45	pWE15, Lawrist	1981	packaging and infection	colonies	5-50
Fosmid	35-45	pFOS1	1992	packaging and infection	colonies	1-2
P1	30-90	pAd10SacBI I	1990	packaging and infection	colonies	1-2
PAC	30-300	pCYPAC2	1994	electroporation	colonies	1-2
BAC	30-300	pBeloBAC1 1	1992	electroporation	colonies	1-2

Cloning Vector: pUC18 Plasmid Vector



- Cloning vector:
1. origin of replication
 2. antibiotics (Amp) resistance gene (=β-lactamase)
 3. lacZ gene (= β-galactosidase)
 4. Multiple cloning site (MCS) at lacZ gene

Recombinant Selection: LacZ & Antibiotics




Blue colony:	non-recombinant
White colony:	recombinant

Sizes of Human Genomic Libraries Prepared in Different Types of Cloning Vector

Type of vector	Insert size (kb)	No. of clones	
		$P = 95\%$	$P = 99\%$
<u>Cosmid</u>	40	240,000	370,000
BAC	300	32,000	50,000
YAC	600	16,000	24,500

(3) 클로닝 기술

Experimental Outline of DNA Cloning






PCR 증폭기술



Polymerase Chain Reaction (PCR)



- Principle of PCR
- Discovery of Taq Enzyme
- Development of Hot start PCR
- Many applications of PCR:
 - RT-PCR
 - RACE
 - DD RT-PCR, etc

Primer Design

예제) Genbank의 p53 유전자 서열

```
1  gtctagagcc accgtccagg gagcaggtag ctgctgggct ccggggacac tttgcgttcg
61  ggctggggagc gtgctttcca cgacggtgac acgcttcctt ggattggcag ccagactgcc
121 ttccgggtca ctgccatgga ggagccgcag tcagatccta gcgtcgagcc ccctctgagt
181 caggaaacat tttcagacct atggaaacta cttcctgaaa acaacgttct gtcccccttg
241 ccgtcccaag caatggatga tttgatgctg tccccggacg atattgaaca atggttcact
301 gaagaccagc gtccagatga agctcccaga atgccagagg ctgctcccc cgtggccctt
361 gcaccagcag ctctacacc ggcgccccct gcaccagccc cctcctggcc cctgtcatct
421 tctgtccctt cccagaaaac ctaccagggc agctacggtt tccgtctggg cttcttgcat
481 tctgggacag ccaagtctgt gacttgcaag tactcccctg ccctcaacaa gatgttttgc
541 caactggcca agacctgccc tgtgcagctg tgggttgatt ccacaccccc gcccggcacc
601 cgcgtccgcg ccatggccat ctacaagcag tcacagcaca tgacggaggt tgtgaggcgc
661 tgccccacc atgagcgctg ctcagatagc gatggtctgg cccctcctca gcattctatc
721 cgagtggaag gaaatttgcg tgtggagtat ttggatgaca gaaacacttt tcgacatagt
781 gtggtggtgc cctatgagcc gcctgagggt ggctctgact gtaccaccat ccactacaac
841 tacatgtgta acagttcctg catgggcggc atgaaccgga ggcccatcct caccatcatc
901 aactggaag actccagtgg taatctactg ggacggaaca gctttgaggt gcgtgtttgt
961 gcctgtcctg ggagagaccg gcgcacagag gaagagaatc tccgcaagaa aggggagcct
1021 caccacgagc tgcccccagg gagcactaag cgagcactgc ccaacaacac cagctcctct
1081 cccagccaa agaagaaacc actggatgga gaatatttca cccttcagat ccgtgggcgt
1141 gagcgcttcg agatgttccg agagctgaat gaggccttgg aactcaagga tgcccaggct
1201 gggaaggagc caggggggag cagggtcac tccagccacc tgaagtccaa aaagggtcag
1261 tctacctccc gccataaaaa actcatgttc aagacagaag ggctgactc agactga
```

Primer Design

- ▶ 길이 20-25mer
- ▶ T_m (melting temp.)
- ▶ GC content
- ▶ Hair pin 구조 유무
- ▶ Self dimer 형성 유

Length

Usually about 20 nt for target sequences in complex genomic DNA; can be much less if target DNA is less complex

Base composition

Substantial tandem repeats of one or more nucleotides to be avoided.

Overall %GC plus length to be chosen so that the T_m of each oligonucleotide (*Table 5.2*) should be equal or nearly identical

Secondary structure

Avoid sequences prone to secondary structure which could form hairpins etc. (see *Figure 1.7A*)

3' end

Base complementarity of the two bases at the extreme 3' end of the two primers to be avoided. Otherwise primer dimers can result, reducing amplification efficiency

<http://www.ncbi.nlm.nih.gov/tools/primer-blast/>

Primer-BLAST

► Primer Design

NCBI/ **Primer-BLAST**: Finding primers specific to your PCR template (using Primer3 and BLAST). [more...](#) [Tips for finding specific primers](#)

PCR Template [Reset page](#) [Save search parameters](#) [Retrieve recent results](#)

Enter accession, gi, or FASTA sequence (A refseq record is preferred) [Clear](#)

Range

From To [Clear](#)

Forward primer

Reverse primer

Use my own forward primer
(5'→3' on plus strand)

Use my own reverse primer
(5'→3' on minus strand)

PCR product size

of primers to return

Primer melting temperatures
(T_m)

Exon/intron selection

Min

Max

70

1000

Min

Opt

Max

Max T_m difference

57.0

60.0

63.0

3

[Add more organisms](#)

Exclusion (optional)

☐ Exclude predicted Refseq transcripts (accession with XM, XR prefix) ☐ Exclude uncultured/environmental sample sequences

Entrez query (optional)

Primer specificity stringency Primer must have at least 2 total mismatches to unintended targets, including

at least 2 mismatches within the last 5 bps at the 3' end.

Ignore targets that have 6 or more mismatches to the primer.

Misprimed product size deviation 4000

Splice variant handling ☐ Allow primer to amplify mRNA splice variants (requires refseq mRNA sequence as PCR template input)

[Get Primers](#)

☐ Show results in a new window ☐ Use new graphic view

Primer 주문

▶ Oligo 주문 form


	Synthesis scale	25 nmole ▼	Purification	HAP ▼
--	-----------------	------------	--------------	-------

● Sequence 입력 * 주문수량 : 개


1	Oligo Name : <input type="text"/>
	Sequence (5'→3') : 5' <input type="text"/>

* 대량주문 의뢰시 페이지 하단의 '대량주문용 양식'을 내려받아 작성후 파일첨부를 이용하시면 편리합니다
* Oligo name, Sequence 입력은 복사후 붙여넣기 하시면 편리합니다.
* Modification Oligo의 경우 추가요청사항에 표기 또는 작성하신 파일을 첨부해 주세요,
* 파일첨부 의뢰시 Oligo Name에 "파일첨부"라고 입력해주세요
* Oligo name 입력시 "W" 는 "WW" 로 표시 됩니다. - "w" 사용시 추가 요청사항에 name을 넣어주세요, -

Workflow: DNA Preparation to Applications




Polymerase Chain Reaction



<http://www.youtube.com/watch?v=2KoLnlwoZKU>

Exponential Amplification of PCR Product



PCR Amplification Protocol

PCR Amplification Protocol


10X PCR buffer containing 15 mM MgCl ₂	2	ul
dNTP (2.5 mM each)	0.8	ul
Forward PCR primer (2.5 uM)	1.2	ul
Reverse PCR primer (2.5 uM)	1.2	ul
template (genomic DNA; 20 ng/4ul)	<u>4</u>	<u>ul</u>
AmpliTaq Gold(5 unit/ul) (1 unit/20 ul Rx)	0.15	ul
dsH ₂ O	10.8	ul

final vol.	20	ul

	94 °C --- 5 min	(1 cycle)
denaturation:	94 °C --- 45 sec	
annealing:	55 °C --- 30 sec	(35 cycles)
extension:	72 °C --- 1:30 sec	
	72 °C --- 10 min	(1 cycle)



PCR Virtual Lab






begin

© 2008 GENETIC SCIENCE LEARNING CENTER, UNIVERSITY OF UTAH

► <http://learn.genetics.utah.edu/content/labs/pcr/>

<http://learn.genetics.utah.edu/content/labs/pcr/>

Microbiome






Types of Cells




▶ Prokaryote 원핵세포

- cells with no nucleus or organelles with membranes.
- Bacteria




▶ Eukaryote 진핵세포


- cells that contain a nucleus and organelles surrounded by a membrane.
- protozoa, algae, fungi, plants, and animals, human




PROCARYOTES





EUCARYOTES





M.T. Madigan and JM Martinko, Brock Biology of Microorganisms を改定









▶ 16S rRNA gene


- ▶ illustrating the conserved (green) and variable



CONSERVED REGIONS: unspecific applications

VARIABLE REGIONS: group or species-specific applications

Bacteria, 유산균... 등



16S-27F : 5' - AGAGTTTGATCCTGGCTCAG-3'

16S-1492R : 5' - GGTTACCTTGTTACGACTT-3'



Fungi, Yeast 동정... 등



ITS1 : 5' - TCCGTAGGTGAACCTGCGG - 3'


ITS4 : 5' - TCCTCCGCTTATTGATATGC - 3'

▶ Universal phylogenetic tree based on the 16S rRNA gene sequence comparisons



칼 우즈 - 1977년 16S 리보솜 RNA의 계통 분류를 통해 새로운 분류군인 고세균을 처음 정의

Jill E. Clarridge III Clin. Microbiol. Rev. 2004;17:840-862



- *E. coli* (*Escherichia coli* DH5a) colony
- White tip을 이용 colony picking 하여
- 10ul의 autoclaved DDW에 희석
- Template로 1ul 사용



Colony PCR




▶ Sample prep


- ▶ Autoclaved DDW 10ul에 picking 한 현탁한 다음 1ul를 template로 사용



▶ EmeraldAmp GT PCR Master Mix

- ▶ Taq polymerase, dNTP, loading dye가 포함되어 있는 mixture 제품.
- ▶ Template, primer 만 더 넣어주면 됨.







○○○ PCR product sequencing 후 ○○○

- ▶ ab1, seq 파일 분석 - FinchTV,
BIOEDIT
- ▶ Alignment
- ▶ Blast Search
- ▶ Phylogenetic Tree





- ▶ 확장자 ab1 파일
- ▶ ABI format






Version: --



Go to Base No. Find Sequence

A C G T

Vertical Scale



Reset Scales

Horizontal Scale



▶ TGCAAGTCGAACGGTAACAGGAAGCAGCTTGCTTCTTTGCTGACGAGTGGCGGACGGGT
GAGTAATGTCTGGGAAACTGCCTGATGGAGGGGGGATAACTACTGGAAACGGTAGCTAAT
ACCGCATAACGTTCGCAAGACCAAAGAGGGGGACCTTCGGGCCTCTTGCCATCGGATGTG
CCCAGATGGGATTAGCTAGTAGGTGGGGTAACGGCTCACCTAGGCGACGATCCCTAGCTG
GTCTGAGAGGATGACCAGCCACACTGGAAGTGAAGACACGGTCCAGACTCCTACGGGAGG
CAGCAGTGGGGAATATTGCACAATGGGCGCAAGCCTGATGCAGCCATGCCGCGTGTATG
AAGAAGGCCCTTCGGGTTGTAAAGTACTTTTCAGCGGGGAGGAAGGGAGTAAAGTTAATAC
CTTTGCTCATTGACGTTACCCGCAGAAGAAGCACCGGCTAACTCCGTGCCAGCAGCCGCG
GTAATACGGAGGGTGCAAGCGTTAATCGGAATTACTGGGCGTAAAGCGCACGCAGGCGG
TTTGTTAAGTCAGATGTGAAATCCCCGGGCTCAACCTGGGAACTGCATCTGATACTGGCA
AGCTTGAGTCTCGTAGAGGGGGGGTAGAATTCAGGTGTAGCGGTGAAATGCGTAGAGAT
CTGGAGGAATACCGGTGGCGAAGGCGGCCCTGGACGAAGACTGACGCTCAGGTGCGA
AAGCGTGGGGAGCAAACAGGATTAGATACCCTGGTAGTCCACGCCGTAAACGATGTCGA
CTTGGAGGTTGTGCCCTTGAGGCGTGGCTTCCGGAGCTAACGCGTTAAGTCGACCGCCTG
GGGAGTACGGCCGCAAGGTTAAAACTCAAaTGAATTGACgGGGGCCCGCACAAAGCGGTG
GAGCATGTGGTTTAATTTCGATGCAACGCGAAGAACCTTACCTGGTCTTGACATCCACGGA
AGTTTTTCAGAGATGAGAATGTGCCTTCGGGAACCGTGAGACAGGTGCTGCATGGCTGTC
GTCAGCTCGTGTTGTGAAATGTTGGGTAAAGTCCCGCAACGAGCGCAACCCTTATCCTTT
GTTGCCAGCGGTCCGGCCGGGAACTCAAAGGAGACTGCCAGTGATAAACTGGAGGAAG
GTGGGGATGACGTCAAGTCATCATGGCCCTTACGACCAGGGCTACACACGTGCTACAATG
GCCCATACAAAGAGAAGCGACCTCGCGAGAGCAAGCGGACCTCATAAAGTGCGTCGTAG
TCCGGATTGGAGTCTGCAACTCGACTCCATGAAGTCGGAATCGCTAGTAATCGTGGATCA
GAATGCCACGGTGAATACGTTCCCGGGCCTTGACACACCGCCCGTCACACCATGGGAGT
GGGTTGCAAAAGAAGTAGGTAGCTTAACCTTCGGGAGGGGCGCTACCAC




BLAST


<http://www.ncbi.nlm.nih.gov/blast>

NCBI/ BLAST/ **blastn**

http://blast.ncbi.nlm.nih.gov/Blast.cgi?PROGRAM=blastn&BLAST_PROGRAMS=megaBlast&PAGE_TYPE=BlastSearch&SHOW_DEFAULTS=on&LINK_LOC=blasthome




SEQ file
BLAST.xlsx
Contig.docx
NCBI Blast(E_coli).htm
treeView(E.coli).png



Ecoli-27F_G09.ab1
Ecoli-27F_G09.seq
Ecoli-1492R_H09.ab1
Ecoli-1492R_H09.seq

- ▶ Seq 파일, 또는
Contig.docx 파일의 염기서열을 이용하여 상동성 분석



	Description	Max score	Total score	Query cover	E value	Ident	Accession
<input type="checkbox"/>	Escherichia coli gene for 16S rRNA, partial sequence, strain: JCM 24009	2604	2604	100%	0.0	99%	AB548582.1
<input type="checkbox"/>	Escherichia coli strain LW1655F+ 16S ribosomal RNA gene, partial sequence	2604	2604	100%	0.0	99%	AY616658.1
<input type="checkbox"/>	Escherichia coli strain HM01 16S ribosomal RNA gene, partial sequence	2603	2603	100%	0.0	99%	JN811622.1
<input type="checkbox"/>	Escherichia coli strain RW-29 16S ribosomal RNA gene, partial sequence	2603	2603	99%	0.0	99%	DQ182324.1
<input type="checkbox"/>	Escherichia coli strain 2 16S ribosomal RNA gene, partial sequence	2601	2601	100%	0.0	99%	JQ907528.1
<input type="checkbox"/>	Escherichia coli strain H1 16S ribosomal RNA gene, partial sequence	2601	2601	100%	0.0	99%	FJ949577.1
<input type="checkbox"/>	Escherichia coli PMV-1 main chromosome, complete genome	2599	18140	100%	0.0	99%	HG428755.1
<input type="checkbox"/>	Escherichia coli LY180, complete genome	2599	18051	100%	0.0	99%	CP006584.1



Copyright © 1997-2013

Tom Hall

Ibis Biosciences

Carlsbad, CA 92008


Biological sequence alignment editor for Win95/98/NT/2K/XP/7
<http://www.ibisbiosciences.edu/BioEdit/bk.html>

- ▶ Contig align
- ▶ BlastN search
- ▶ Genetic Tree



Reverse Transcriptase (RT) PCR 기술

Workflow of RT-PCR: RT Reaction for cDNA Synthesis from Total RNA



Real-time or Quantitative PCR (qPCR): Measuring Gene Expression Levels

(A) TaqMan Assay


Annealing

R is reporter fluorophore, which emits at a wavelength absorbed by the quencher fluorophore (Q).




Probe displacement

DNA polymerase starts extending primers moving toward the probe.




Probe cleavage

The probe is degraded. The reporter is released from the quencher and starts to emit fluorescence.

(B) Relative Quantitation: $C_t - \Delta C_t - \Delta \Delta C_t$ 

C_t = the cycle No. at which the fluorescence signal crosses the threshold.



$$\Delta C_{t0} = \text{target} - \text{ref} = 9.70$$

$$\begin{aligned} \text{Difference } (\Delta \Delta C_t) &= \Delta C_{t1} - \Delta C_{t0} \\ &= 9.70 - (-1.7) = 11.40 \end{aligned}$$


$$\begin{aligned} \text{RQ (relative quantitation)} &= \text{fold change} = 2^{-\Delta \Delta C_t} \\ &= 2^{11.40} = 2702 \end{aligned}$$

$$\Delta C_{t1} = \text{target} - \text{ref} = -1.7$$

- C_t is inversely correlated with starting copies (ie greater C_T = less starting copies).
- C_T depends on where the threshold is set.
- Same threshold must be set across samples for meaningful data comparison.

PCR vs Real-time PCR

► Quantitative real-time PCR



Principle of Digital PCR

Digital PCR (dPCR) is a refinement of conventional PCR methods that can be used to directly quantify and clonally amplify nucleic acids (including DNA, cDNA, methylated DNA, or RNA). The key difference between dPCR and traditional PCR lies in the method of measuring nucleic acids amounts, with the former being a more precise method than PCR. PCR carries out one reaction per single sample. dPCR also carries out a single reaction within a sample, however the sample is separated into a large number of partitions and the reaction is carried out in each partition individually. This separation allows a more reliable collection and sensitive measurement of nucleic acid amounts. The method has been demonstrated as useful for studying variations in gene sequences - such as copy number variants, point mutations, and it is routinely used for clonal amplification of samples for "next-generation sequencing."

Conventional PCR



Split sample by dilution

Digital PCR



— wild type
— mutant



One measurement

Droplet digital PCR




Many thousands of
many measurements

- **RainDance** (1,000,000 droplets)
- **BioRad** (20,000 droplets)
- **Life Tech** (20,000 chambers)
- **Fluidigm**

Applications of Digital PCR

- **Copy number variation (CNV)** — 1.2-fold ($\pm 10\%$) differences in gene copy numbers
- **Rare sequence detection** — as low as 1 in 100,000 (for example, a few tumor cells in a wild-type background)
- **Gene expression and miRNA analysis** — absolute quantification of expression levels, especially low-abundance miRNAs.
- **Single cell analysis** — the high degree (10–100-fold) of cell-cell variation in gene expression levels among homogeneous single-cell populations and stem cells drives the need for the analysis of single cells; digital PCR enables low copy number target quantification
- **Pathogen detection** — detection and quantification of pathogens in samples quickly and precisely in very small amounts in biological samples.
- **Next-generation sequencing (NGS)** — quantification of NGS sample preparations to increase sequencing accuracy and reduce run repeats, and validation of sequencing results with absolute quantification


1st



PCR
Qualitative




2nd



Real-time PCR
Relative Quantification



3rd



Digital PCR
Absolute Quantification

